




Product Datasheet

ACF antibody (orb28588)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	A1CF
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65202
Uniprot ID	Q9NQ94
NCBI	NP_001185748.1 , NP_620311.1 , NP_055391.2 , NP_620310.1 , NP_001185749.1 , NP_001185747.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis
of A375 cell line ...



Immunohistochemical
analysis of formalin...