

Product Datasheet

Metallothionein 3 antibody (orb2871)

Description	Metallothionein 3 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Metallothionein 3 (1-68/68 aa)
Target	Metallothionein 3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	7.5 kDa
Uniprot ID	P25713
Dilution Range	IHC-P=1:100-500, Flow-Cyt=0.2ug/test, WB=1:500-2000, IHC-F=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.