

Product Datasheet

RPN13 antibody (Biotin) (orb28748)

Description	RPN13 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Saccharomyces cerevisiae 26S proteasome regulatory subunit RPN13 protein (2-156AA)
Target	RPN13
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O13563
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.