



Product Datasheet

VCX1 antibody (orb28760)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	VCX
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	22289
Uniprot ID	Q9H320
NCBI	NP_038480.2
Dilution Range	IHC-P: 1:50-100, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of A2058, MCF-7, M...