

Product Datasheet

IspL antibody (Biotin) (orb28772)

Description	IspL antibody (Biotin)
Species/Host	Rabbit
Reactivity	Bacteria
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Rhizobium meliloti UDP-glucuronate 5\'-epimerase protein (1-341AA)
Target	IspL
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	O54067
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.