




Product Datasheet

UBTD1 antibody (orb28775)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	UBTD1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	25938
Uniprot ID	Q9HAC8
NCBI	NP_079230.1
Dilution Range	WB: 1:1000, IHC-P: 1:50
Expiration Date	12 months from date of receipt.



Western blot analysis of HL-60 cell line...



Immunohistochemical analysis of formalin...