

Product Datasheet

ENT1 antibody (orb28810)

Description Peptide affinity purified rabbit polyclonal antibody.

Species/Host Rabbit

Reactivity Human, Mouse

Conjugation Unconjugated

Tested Applications WB

Immunogen Synthetic Peptide

Target SLC29A1

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal

MW 50219

Uniprot ID [Q99808](#)

NCBI [NP_001071644.1](#), [NP_001071642.1](#),
[NP_001071643.1](#), [NP_004946.1](#), [NP_001071645.1](#)

Dilution Range WB: 1:1000

Expiration Date 12 months from date of receipt.

