



Product Datasheet

AMN1 antibody (orb28844)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	AMN1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	28408
Uniprot ID	Q8IY45
NCBI	NP_001265340.1 , NP_001106873.1 , NP_001265341.1
Dilution Range	FACS: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of K562, U251, A54...