




Product Datasheet

CSH1 antibody (orb28872)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	This CSH1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 177-204 amino acids from the C-terminal region of human CSH1.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Uniprot ID	P01243
NCBI	NP_001308.1 , NP_066271.1
Dilution Range	WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cell lysate...



Immunohistochemical analysis of formalin...