

Product Datasheet

KIF5A antibody (orb2888)

Description	KIF5A antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human KIF5A/NKHC1 (201-300/1032 aa)
Target	KIF5A/NKHC1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	113 kDa
Uniprot ID	Q12840
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.