

Product Datasheet

pCPLA2 antibody (orb2891)

Description	pCPLA2 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human CPLA2 around the phosphorylation site of Ser505 (PL(p-S)PL)
Target	phospho-CPLA2 (Ser505)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	82 kDa
Uniprot ID	P47712
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.