



## Product Datasheet

### NUSAP antibody (orb28942)

|                            |                                                                                                                                                                              |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Peptide affinity purified rabbit polyclonal antibody.                                                                                                                        |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                       |
| <b>Reactivity</b>          | Human                                                                                                                                                                        |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                 |
| <b>Tested Applications</b> | IHC-P, WB                                                                                                                                                                    |
| <b>Immunogen</b>           | Synthetic Peptide                                                                                                                                                            |
| <b>Target</b>              | NUSAP1                                                                                                                                                                       |
| <b>Form/Appearance</b>     | Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification. |
| <b>Note</b>                | For research use only                                                                                                                                                        |
| <b>Isotype</b>             | Rabbit IgG                                                                                                                                                                   |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                   |
| <b>MW</b>                  | 49452                                                                                                                                                                        |
| <b>Uniprot ID</b>          | <a href="#">Q9BXS6</a>                                                                                                                                                       |
| <b>NCBI</b>                | <a href="#">NP_001230072.1</a> , <a href="#">NP_001230071.1</a> ,<br><a href="#">NP_057443.2</a> , <a href="#">NP_001230073.1</a> , <a href="#">NP_060924.4</a>              |
| <b>Dilution Range</b>      | WB: 1:1000, IHC-P: 1:50-100                                                                                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                              |



Immunohistochemical analysis of formalin...