




Product Datasheet

HDAC2 antibody (orb28974)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	HDAC2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	55364
Uniprot ID	Q92769
NCBI	NP_001518.3
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



WB analysis of 293 cell line lysate using...



WB analysis of mouse spleen tissue lysate...