

Product Datasheet

BCDO2 antibody (orb28999)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	BCO2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	65674
Uniprot ID	Q9BYV7
NCBI	NP_001243327.1 , NP_001243326.1 , NP_001032367.2 , NP_001243329.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

