



Product Datasheet

UNCX antibody (orb29129)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	UNCX
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	53690
Uniprot ID	A6NJT0
NCBI	NP_001073930.1
Dilution Range	IF/ICC: 1:25, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of U251 cell...