



Product Datasheet

RNF19B antibody (orb29145)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	RNF19B
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	77925
Uniprot ID	Q6ZMZ0
NCBI	NP_001120833.1 , NP_699172.2
Dilution Range	WB: 1:500, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of RNF19B antibody...