



Product Datasheet

LARS antibody (orb29172)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	LARS1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	134466
Uniprot ID	Q9P2J5
NCBI	NP_064502.9
Dilution Range	FACS: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of LARS antibody i...