

Product Datasheet

SMT3 antibody (Biotin) (orb29250)

Description	SMT3 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Saccharomyces cerevisiae Ubiquitin-like protein SMT3 protein (2-98AA)
Target	SMT3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q12306
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.