

## Product Datasheet

### PDE6H antibody (FITC) (orb29266)

|                        |                                                                                                                       |
|------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | PDE6H antibody (FITC)                                                                                                 |
| <b>Species/Host</b>    | Rabbit                                                                                                                |
| <b>Reactivity</b>      | Human                                                                                                                 |
| <b>Conjugation</b>     | FITC                                                                                                                  |
| <b>Immunogen</b>       | Recombinant Human Retinal cone rhodopsin-sensitive cGMP 3',5'-cyclic phosphodiesterase subunit gamma protein (1-83AA) |
| <b>Target</b>          | PDE6H                                                                                                                 |
| <b>Preservatives</b>   | Preservative: 0.03% Proclin 300.<br>Constituents: 50% Glycerol                                                        |
| <b>Form/Appearance</b> | Liquid                                                                                                                |
| <b>Storage</b>         | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                         |
| <b>Note</b>            | For research use only                                                                                                 |
| <b>Isotype</b>         | IgG                                                                                                                   |
| <b>Clonality</b>       | Polyclonal                                                                                                            |
| <b>Uniprot ID</b>      | <a href="#">Q13956</a>                                                                                                |
| <b>Dilution Range</b>  | ELISA                                                                                                                 |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                       |