

Product Datasheet

AP1M1 antibody (orb29319)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Synthetic Peptide
Target	AP1M1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	48587
Uniprot ID	Q9BXS5
NCBI	NP_115882.1 , NP_001123996.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

