


Product Datasheet

Collagen VI antibody (orb2936)

Description	Collagen VI antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human COL6A5/Collagen VI alpha 5 (121-220/2615 aa)
Target	COL6A5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	288 kDa
Uniprot ID	A8TX70
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...