




Product Datasheet

ATG4D antibody (orb29416)

Description	Purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	ATG4D
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	52922
Uniprot ID	Q86TL0
NCBI	NP_001268433.1 , NP_116274.3
Dilution Range	IF/ICC: 1:200, WB: 1:1000, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse liver tis...



Immunohistochemical analysis of formalin...