



Product Datasheet

C1QL4 antibody (orb29483)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	C1QL4
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	24909
Uniprot ID	Q86Z23
NCBI	NP_001008224.1
Dilution Range	FACS: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
C1QL4
antibo...