

Product Datasheet

TUBA3C antibody (HRP) (orb29539)

Description	TUBA3C antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Tubulin alpha-3C/D chain protein (161-418AA)
Target	TUBA3C
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13748
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.