




Product Datasheet

Osteopontin antibody (orb29541)

Description	Rabbit polyclonal antibody to Osteopontin
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	SPP1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35423
Uniprot ID	P10451
NCBI	NP_001035147.1, NP_001238759.1, NP_001035149.1, NP_000573.1, NP_001238758.1
Dilution Range	WB: 1:1000, IF/ICC: 1:25, IHC-P: 1:50-100
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...