

Product Datasheet

LAPTM5 antibody (Biotin) (orb29606)

Description	LAPTM5 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Lysosomal-associated transmembrane protein 5 protein (205-262AA)
Target	LAPTM5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q13571
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.