

Product Datasheet

LILRA5 antibody (orb29617)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LILRA5
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	32755
Uniprot ID	A6NI73
NCBI	NP_870994.1 , NP_871714.1 , NP_067073.1 , NP_871715.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

