




Product Datasheet

HRAS antibody (orb29649)

Description	Synthetic peptide.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Synthetic Peptide
Target	HRAS
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is prepared by Saturated Ammonium Sulfate (SAS) precipitation followed by dialysis against PBS.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	21298
Uniprot ID	P01112
NCBI	NP_001123914.1 , NP_789765.1 , NP_005334.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse cerebellu...



Confocal immunofluorescence analysis of ...