



Product Datasheet

OC90 antibody (orb29651)

Description	Rabbit polyclonal antibody to LSM7 is a member of the LSM subfamily of snRNP Sm proteins. LSM7 is in...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	OC90
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	51728
Uniprot ID	Q02509
NCBI	NP_001073868.2
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
K562 cells (...)