

Product Datasheet

LSM7 antibody (orb29652)

Description	Rabbit polyclonal antibody to NCL is a major nucleolar phosphoprotein that is involved in pre-RNA tr...
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	LSM7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	11602
Uniprot ID	Q9UK45
NCBI	NP_057283.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

