



Product Datasheet

CPA3 antibody (orb29660)

Description	Rabbit polyclonal antibody to PALB2 is significant to the BRCA2 function, DNA damage response and ch...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	Synthetic Peptide
Target	CPA3
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	48670
Uniprot ID	P15088
NCBI	NP_001861.2
Dilution Range	WB: 1:1000, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow
cytometric
analysis of
K562 cells (...)