




Product Datasheet

QTRTD1 antibody (orb29669)

Description	Rabbit polyclonal antibody to SSH3 is localized at both the nucleus and cytoplasm. SSH2 performs cat...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	QTRT2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	46713
Uniprot ID	Q9H974
NCBI	NP_001243765.1 , NP_001243766.1 , NP_001243764.1 , NP_078914.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50, IF/ICC: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Confocal immunofluorescence analysis of ...