




Product Datasheet

CREB1 antibody (orb29680)

Description	Rabbit polyclonal antibody to PTPRC belongs to the protein-tyrosine phosphatase family. This is a re...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, IF, WB
Immunogen	Synthetic Peptide
Target	CREB1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	35136
Uniprot ID	P16220
NCBI	NP_604391.1 , NP_004370.1
Dilution Range	WB: 1:1000, IF/ICC: 1:10-50, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of A549 cells (...)



Confocal immunofluorescence analysis of ...