




Product Datasheet

TMEM100 antibody (orb29686)

Description	Rabbit polyclonal antibody to BRUNOL6 is a member of the CELF/BRUNOL family. It is expressed primari...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-P, WB
Immunogen	Synthetic Peptide
Target	TMEM100
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	14386
Uniprot ID	Q9NV29
NCBI	NP_001093110.1 , NP_060756.2
Dilution Range	WB: 1:1000, IHC-P: 1:50-100, FACS: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of NCI-H460 cel...