

Product Datasheet

TTC6 antibody (orb29697)

Description	Synthetic peptide.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TTC6
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	59379
Uniprot ID	Q86TZ1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

