



Product Datasheet

KRCC1 antibody (orb2972)


Description	KRCC1 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human KRCC1 (161-259/259 aa)
Target	KRCC1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	31 kDa
Uniprot ID	Q9NPI7
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of mouse kidney tissue.



Western blot analysis of mouse kidney lysates.



Immunohistochemical staining of Mouse brain tissue.