



Product Datasheet

CCDC110 antibody (orb29730)

Description	Rabbit polyclonal antibody to PALM3 is a member of the paralemmin family that localizes to the cyto...
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	CCDC110
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	96726
Uniprot ID	Q8TBZ0
NCBI	NP_001138883.1 , NP_689988.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...