

Product Datasheet

PALM3 antibody (orb29731)

Description Rabbit polyclonal antibody to TAS1R2 is a member of the G-protein coupled receptor 3 family. It is ...

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications FC, IHC-P, WB

Immunogen Synthetic Peptide

Target PALM3

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 71695

Uniprot ID [A6NDB9](#)


NCBI [NP_001138500.1](#)

Dilution Range WB: 1:1000, IHC-P: 1:10-50, FACS: 1:10-50

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of formalin...



Flow cytometric analysis of A549 cells (...)