

Product Datasheet

TAS1R2 antibody (orb29732)

Description	Rabbit polyclonal antibody to LRP12, which is also known as ST7 belongs to LDLR family. It is invol...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	TAS1R2
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	95183
Uniprot ID	Q8TE23
NCBI	NP_689418.2
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

