



Product Datasheet

RHOJ antibody (orb29735)

Description	Rabbit polyclonal antibody to SET7 belongs to the SET domain-containing family that is highly activ...
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	RHOJ
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	23821
Uniprot ID	Q9H4E5
NCBI	NP_065714.1
Dilution Range	IHC-P: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...