


Product Datasheet

KIF1B antibody (orb2974)

Description	KIF1B antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human KIF1B (13-120/1816 aa)
Target	KIF1B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	204 kDa
Uniprot ID	O60333
Expiration Date	12 months from date of receipt.



Paraffin-
embedded rat
brain tissue
fixed...