

Product Datasheet

NUP107 antibody (orb29759)

Description Rabbit polyclonal antibody to NUP107 also known as Nucleoporin Nup107 whose expression can be ubiqui...

Species/Host Rabbit

Reactivity Mouse

Conjugation Unconjugated

Tested Applications IHC-P, WB

Immunogen Synthetic Peptide

Target NUP107

Form/Appearance Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.

Note For research use only

Isotype Rabbit IgG

Clonality Polyclonal



MW 106374

Uniprot ID [P57740](#)

NCBI [NP_065134.1](#)

Dilution Range WB: 1:1000, IHC-P: 1:50-100

Expiration Date 12 months from date of receipt.



Immunohistochemical analysis of formalin...