




Product Datasheet

DEF antibody (orb29782)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Hamster, Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	UTP25
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	87055
Uniprot ID	Q68CQ4
NCBI	NP_055203.4
Dilution Range	WB: 1:1000, IHC-P: 1:10-50, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...



Western blot analysis of CHO cell line I...