

Product Datasheet

tcag7.1227 antibody (orb29813)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	This tcag7.1227 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 183-211 amino acids from the Central region of human tcag7.1227.
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.

