




Product Datasheet

PRKAG1 antibody (orb29825)

Description	Rabbit polyclonal antibody to PRKAG1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Target	PRKAG1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	37579
Uniprot ID	P54619
NCBI	NP_002724.1 , NP_001193639.1 , NP_001193638.1
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse bladder t...



Western blot analysis of 293 cell lysate...