


Product Datasheet

ZBT24 antibody (orb29856)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	ZBTB24
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	78282
Uniprot ID	O43167
NCBI	NP_055612.2 , NP_001157785.1
Dilution Range	WB: 1:1000, IHC-P: 1:10-50
Expiration Date	12 months from date of receipt.



Western blot analysis of ZBT24 antibody ...



Western blot analysis of ZBT24 antibody ...