

Product Datasheet

HFE antibody (orb29868)

Description	HFE antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human Hereditary hemochromatosis protein (23-306AA)
Target	HFE
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q30201
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.