




Product Datasheet

WT1 antibody (orb29938)

Description	Rabbit polyclonal antibody to WT1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	FC, IF, IHC-P, WB
Immunogen	Synthetic Peptide
Target	WT1
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	49188
Uniprot ID	P19544
NCBI	NP_001185481.1 , NP_000369.3 , NP_077742.2 , NP_001185480.1 , NP_077744.3
Dilution Range	FACS: 1:10-50, WB: 1:1000, IHC-P: 1:50-100, IF/ICC: 1:10-50
Expiration Date	12 months from date of receipt.



Confocal immunofluorescent analysis of M...



Western blot analysis of WT1 antibody in...