



Product Datasheet

TRIM7 antibody (orb29952)

Description	Peptide affinity purified rabbit polyclonal antibody.
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	Synthetic Peptide
Target	TRIM7
Form/Appearance	Purified polyclonal antibody supplied in PBS with 0.09% (W/V) sodium azide. This antibody is purified through a protein A column, followed by peptide affinity purification.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56631
Uniprot ID	Q9C029
NCBI	NP_976039.1 , NP_976038.1 , NP_976042.1 , NP_976041.1 , NP_203128.1 , NP_976040.1
Dilution Range	IHC-P: 1:10-50, WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of NCI-H292