

Product Datasheet

KRTAP19-5 antibody (HRP) (orb29953)

Description	KRTAP19-5 antibody (HRP)
Species/Host	Rabbit
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Keratin-associated protein 19-5 protein (1-72AA)
Target	KRTAP19-5
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q3LI72
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.